

INTRODUCTION


We are creating a database of pulse waveforms for scintillation detectors over gamma-ray energies from 30-1000 keV. This work will aid research teams around the globe who are focused on building the next generation of instruments for gamma-ray astronomy.

This poster provides initial results from the database using a Hamamatsu R12699 flat panel photomultiplier tube (PMT) coupled to NaI(Tl) and CsI(Na) scintillators.

BACKGROUND & SETUP


- ❖ A scintillator crystal emits photons in the visible spectrum when hit with high-energy photons [1].
- ❖ A photosensor accepts these photons and sends an analog signal to the digitizer.
 - Our photosensor is a Hamamatsu R12699 flat panel PMT.
- ❖ The digitizer then outputs a digital signal of the waveform to the computer to be saved.

Figure 1: Lab Setup for the Experiment




METHODS

- ❖ For our work, we use a “phoswich” style detector, which couples two scintillators to a single photosensor, because it allows us to acquire waveforms for two scintillators at once.
- ❖ The first graph returned (Figure 2.1) is an energy histogram of Ba-133 for both scintillators in the NaI(Tl)/CsI(Na) phoswich, which we then need to separate based on individual pulse shapes.



- ❖ We separate energy values for each scintillator by integrating the pulse shape over two timescales, defined by long and short energy gates. Following reference [2], we then formulate a pulse shape discrimination (PSD) variable as a ratio and produce the graph displayed in Figure 2.2.





$$PSD = \frac{\text{Long Gate} - \text{Short Gate}}{\text{Long Gate}}$$

- Smaller values of PSD correspond to faster scintillator pulse shapes.

- ❖ Selecting energy values from the population of smaller PSD values yields a spectrum for NaI(Tl) alone. We can see that there are three major peaks in the light curve (Figure 2.3). These peaks correspond to Ba-133's most common energy levels [3]:


- 30 keV
- 81 keV
- 356 keV

- ❖ Sorted waveform data is used to plot the average waveform at the energy parameters for the desired scintillator (Figure 2.4).



RESULTS & CONCLUSIONS

- ❖ We can compare the waveforms of the different scintillators at the same energy to estimate the minimum energy thresholds (Figure 3).
- ❖ For instance, since CsI(Na) has a lower peak amplitude than NaI(Tl) at the same energy, only NaI(Tl) would register for a threshold set at the level of the blue line.
 - So, if we wanted to capture energy events higher than 0.05 V, a NaI(Tl) scintillator would be desired over CsI(Na).



CITATIONS

- [1] G. Knoll, *Radiation Detection and Measurement*, 3rd ed., 2000.
- [2] UM5960 - CoMPASS User Manual, Revision 14, CAEN, 2021.
- [3] S.Y.F. Chu et al., *The Lund/LBNL Nuclear Data Search*, ver. 2.0, 1999, <http://nucleardata.nuclear.lu.se/toi/perchart.htm>.

ACKNOWLEDGEMENTS

This work was supported by the NSF award AGS-1950831. Any opinions, findings, and conclusions or recommendations expressed in this material are those of the author(s) and do not necessarily reflect the views of the National Science Foundation.

I'd also like to thank Dr. Wood for his support in this project.